Bilancia Digitale

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Package contents

Scale
Alkaline batteries 2 x 1.5V, type AA
Instruction manual

With its ergonomic and compact shape, the Chicco Digital Scale weighs baby from the first few days with a capacity ranging from 30 g to 20 kg. The 10 g resolution also allows for the measurement of small variations in weight; a sophisticated anti-vibration program stabilises the weight to quickly and accurately measure the weight of baby even when he is moving.

Legend (Fig. A)

1. Scale weighing plate
2. LCD display
3. On/off button 
4. Tare function key 
5. Battery Compartment

WARNINGS

- This scale is only to be used for domestic purposes and may not be used for public use or for commercial transactions.
- Always lift the scale when moving it. Do not drag or scrape it as it may become damaged or its functioning may be impaired.
- Always use the scale with the utmost care.
- Never leave baby on the scale unattended!
- Only rest the scale on flat, stable surfaces.
- Never place weights on the scale that exceed its maximum capacity (20 kg).
- If the scale is not used for long periods, remove the batteries from the battery compartment.
- Keep the components of the packaging out of the reach of children.
- After having removed the appliance from its packaging, ensure that the scale is intact and that there is no visible damage. If in doubt, do not use it and only contact the retailer.
- Do not immerse the scale in water or any other liquid, do not place it under a jet of water or allow water to get into it, or use it in the proximity of sources or accumulations of water: the scale could be damaged or its functioning could be impaired.
- Product features may be subject to change without notice.

HOW TO USE

1. PREPARATION

Remove the battery compartment cover on the bottom of the scale (Fig. 1), insert the batteries supplied, taking care to observe the correct polarity as indicated on the inside of the compartment, and replace the cover correctly. Position the scale on a horizontal smooth and stable surface without dragging it or scraping it along the surface (Fig. 2). Never place weights on the scale that exceed its maximum capacity (20 kg).

2. SWITCHING ON AND DISPLAY OF LAST WEIGHT

Press the on/off key , the display will indicate the symbols  for about 1 second. Then the value of the last recorded weight out will appear for about 2 seconds followed by the indication  (Fig.3). It is possible, thanks to the memory storage of the last weight, to calculate the difference with the weighing in progress: simply make a note of the weight that will appear before performing the next weighing.

2.2 WEIGHT SET CHANGE

After the device is powered ON and when  appears, keep the  button pressed for 5sec in order to switch from KG to LB/OZ (Fig. 4).

2.3 WEIGHT MEASUREMENT

When the scale display indicates , it is possible to place a towel on the weighing plate and then press the TARE key  (Fig. 4). Check that the display indicates , then place baby in the centre of the scale weighing plate (Fig. 5).

ATTENTION: the scale is equipped with an automatic device to stabilise the weight to eliminate the effects of measuring interference caused by the inevitable movements of baby. The waiting time may vary depending on the movement of baby. When the weight has stabilised, the display will flash 3 times, and after the third flash the weight of the baby will appear (Fig. 6). The scale will automatically switch off after 30 seconds (Fig. 7).

ATTENTION: with weights greater than 20 kg, the scale will show the following error message  (Fig. 8).

3. FLAT BATTERIES

The display of the symbols  , indicates the flat status of the batteries which should be replaced.

3.1 REMOVAL, REPLACEMENT AND INSERTION OF THE BATTERIES

- Remove the battery compartment cover on the bottom of the scale (Fig. 9) and remove the batteries. Insert the new batteries ensuring that they are of the same type or equivalent to the ones indicated, making sure to observe the correct polarity as indicated on the inside of the compartment; replace the cover correctly.
- The batteries must only be replaced by an adult.
 - Use alkaline batteries equal or equivalent to those indicated in this manual.
 - Do not mix different types of batteries or flat batteries with new ones.
 - Do not leave batteries or any tools within the reach of children.
 - Do not short-circuit the battery poles.
 - if the product is not used for a long period, always remove the batteries.
 - Always remove flat batteries from the product to avoid any loss of fluid that can damage it.
 - Dispose of flat batteries by recycling, do not throw away, do not burn or dispose of with municipal waste.
 - If the batteries are leaking liquid, replace them immediately, making sure to clean the battery compartment; thoroughly wash your hands if they have come into contact with the liquid.
 - Do not attempt to recharge non-rechargeable batteries: they may explode.
 - Using rechargeable batteries is not advisable as doing so may impair the functionality of the product.
 - The product is not designed to run on replaceable lithium batteries.

WARNING: improper use may create hazardous conditions.

4. CLEANING AND MAINTENANCE

Clean the scale using a damp cloth, do not use abrasive products or chemicals. Do not soak or immerse the scale in water. When not in use, store the scale in a cool dry place away from light and/or heat sources, avoiding placing objects on top of it.

TECHNICAL SPECIFICATIONS

Mod. 06746/A
REF. 00 005577 000 000
Size: mm 610 x 345 x 64 mm
Min. capacity 0.03 kg (30 g)
Max. capacity 20 kg (20,000 g)
Graduation: 10 g
Power supply: 2 x 1.5 V AA alkaline batteries

Operating temperature range¹ +10/ +40°C
Transport and storage temperature range -10/+50°C

FUNCTIONS

Weight lock (Vibration dampers)
Automatic shut-off
Tare button 
Last recorded weight
Low battery indicator
Overload indicator
LCD display

SYMBOLS LEGEND

 Product meets the requirements of the applicable EC directives

 WARNING!

 Read the instructions manual.

 Instructions for use

DECLARATION OF CONFORMITY.

Artsana S.p.A. hereby declares that this Chicco brand scale model 06746/A complies with all essential requirements and other provisions indicated in the EU Directives 2004/108/EC (electromagnetic compatibility) and 2011/65/EU (RoHS). Please feel free to contact Artsana S.p.A. Via S. Catelli 1 – 22070 Grandate (Co) Italy for a copy of this EC Declaration.

 **This product complies with the Directive 2012/19/EU.**
The crossed bin symbol on appliance indicates that the product, at the end of its life, must be disposed of separately from domestic waste, either by taking it to a separate waste disposal site for electric and electronic appliances or by returning it to your dealer when you buy another similar appliance. The user is responsible for taking the appliance to a special waste disposal site at the end of its life. If the disused appliance is collected correctly as separate waste, it can be recycled, treated and disposed of ecologically; this avoids a negative impact on both the environment and health, and contributes towards the recycling of the product's materials. For further information regarding the waste disposal services available, contact your local waste disposal agency or the shop where you bought the appliance.

 **THIS PRODUCT COMPLIES WITH EU DIRECTIVE 2006/66/EC as amended.**
The crossed bin symbol on the batteries indicates that, at the end of their life, they must be disposed of separately from domestic waste, either by taking them to a separate waste disposal site for batteries or by returning them to your dealer when you buy similar rechargeable or non-rechargeable batteries. The chemical symbols Hg, Cd, Pb, printed under the crossed bin symbol, indicate the type of substance contained in the batteries: Hg=Mercury, Cd=Cadmium, Pb=Lead. The user is responsible for taking the batteries to a special waste disposal site at the end of their life so that they can be treated and recycled. If the spent batteries are correctly collected as separate waste, they can be recycled, treated and disposed of ecologically; this avoids a negative impact on both the environment and human health, and contributes towards the recycling of the batteries' substances. Non compliance with battery disposal regulations damages the environment and human health. For further information regarding the waste disposal services available, contact your local waste disposal agency or the shop where you bought the batteries.

ARTSANA reserves the right to make changes to this instruction manual at any time and without warning. The reproduction, transmission, transcription or translation into other languages of this manual, or part of the same, is strictly prohibited without prior written authorisation from ARTSANA.

WARRANTY

The product is guaranteed against any conformity defect in normal conditions of use as provided for by the instructions for use. The warranty shall not therefore apply in the case of damages caused by improper use, wear or accidental events. For the duration of warranty on conformity defects please refer to the specific provisions of national laws applicable in the country of purchase, where provided.

Artsana S.p.A.
Via Saldarini Catelli 1 - 22070 Grandate (Como) Italy
www.chicco.com

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