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Thanks for you to select our product. Before use, please carefully read the manual and pay more attention to safety notice item in order to correctly operate and use.

Please take good care of the manual in convenience that you can query.

SAFETY ITEMS

Use Condition

Please don't use in wet or dirty condition to prevent any trouble from happening. Please don't use this product in room temperature above than 40°C.

Please don't put the product outside room.

Unfit people to use this product

If you have the following state, before use, please consult doctor. Curer who is receiving the cure from doctor or who is feeling bad.

Malignancy Patient

Pregnant Woman or Woman in menses Urgent illness Patient

Controller Notice Items

- Don't use your wet hand to operate the controller.
- Don't fall the controller to the floor or in water.
- Don't put any heavy goods on controller.

NOTICE

During massage, the massage head should be on both sides of the neck. Otherwise maybe it will result in ache.

Please don't sleep when massaging.

Please don't use this produce when drunken or bad.

Please don't directly keep massage wheel near to hand, elbow parts, knee, breast or belly.

Please don't exceed 30 minutes once when continue to use the massage.

DISUSE PRODUCE CONDITION

When there is the following, please Immediately stop using the product Close Controller Switch and Plug out or smoke socket.

When water or other liquid is spilled on the product or the remote control.

Some of the internal parts will be damaged or will be visible on the outside.

During massaging feel abnormal or have any severe ache.

During massage find any abnormal state such as fume, peculiar smell and etc.

Suddenly power broken

Thunder occur

SAFETY ITEMS

Product Assemble &Maintenance.

Please don't privately disassembly the packing of this product.

Do not touch the internal parts of this product to avoid any problems.

This products can be maintained only by professional service worker.

Do not disassemble the device yourself. .

ABOUT POWER ADAPTER &POWER CABLE

Please check whether power adapter accord with the specified voltage, current and match with the output

When unplugging the power cord, grasp the plug, not the cord.

It is strictly forbidden to pull the cord.

Be careful that heavy objects do not pinch the cord.

INNER PART NAME & FUNCTIONS

Appearance Construction

PRODUCT TECHNICAL PARAMETER

Product Name:	Multi-purpose massage pillow
DC Input Voltage:	12 V
Aln put Current:	1. 5- 2 .0 Rating Power: 20W
Safe Design:	[Q]
Power Adapter:	- - lit - +

USE METHOD

Plug controller into the socket of smoke source(In car use) or plug the power adapter in the power socket.

Dial direction switch key, massage pillow enter into massaging state (Advise: Massaging 10-15 minutes /once).

Turn off the device by switching to the middle mode, then disconnect from the socket.

FAQ & SOLUTION

Sequence	Phenomenon	Possible Reason	Solution
1	In operating appear „xiuxiu“ sound. Product in working order.	Or the frictional sound between massage ball and fabric packing	No need any disposal
2	Open to fail	Product is too hot to protect	Close the product for 30 minutes to use it again
		Power Plug or Smoke source socket is not plugged well	Plug the socket well
		Direction Switch dial to middle position	Direction toward left or right

If dispose the above, still there is any abnormal situation please contact with manufacturer and request them to maintain and repair.

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