# CHA-CHA

## DIGITAL FAT SCALE

**MANUAL**

**EBS022K EBS022W**

# DIGITAL FAT SCALE

## Product features

- . It can be used as a normal scale
- . Body fat scale can measure weight, but it is specially designed for measurement of the fat, water and muscle ranges of your body.

## Technical specification

- Measuring range: 5 kg-180 kg/ Weighing sensitivity: 0.1 kg
- Fat measurement range: 5-50% / Fat measurement precision: 0.1%
- Water measurement range: 20-70% / Water measurement precision: 0.1%
- Muscle measurement range: 20-56% / Muscle measurement precision: 0.1%
- Applicable height range: 75-225 cm
- Applicable scope of age: 10-85 years old
- Minimum weighing: 5 kg
- Three units: kg/st/lb
- It can store data for 12 users
- The error indicator: Err %
- Voltage display: the battery voltage is less than 2.5 V, shows "Lo" and 2 seconds later the scale will auto power-off.
- Overload display: O-Ld, overload value: 182 kg
- The auto power-off: No operation about 10 seconds will automatically switch off
- Function buttons: UP, SET, DOWN, UNIT

## Measuring fat/water/muscle mode

### 1 Setting of personal parameters. The parameters may be set in the following ranges:

- Press Key **[SET]** to turn on, the user serial number P00 (P00-P11) flickers, press key **[▲]** or **[▼]** to choose the number which you want (12 numbers in total), then press **[SET]** to confirm. If there is no operation within 3 seconds, system will confirm automatically that the current flickering number is the current user data.
- After user number setting, gender sign starts flickering. Press **[▲]** or **[▼]** to choose gender and press **[SET]** to confirm. If there is no operation, it will be turn off automatically within 20 seconds.
- After gender setting, age sign will start flickering. Press **[▲]** or **[▼]** to adjust age, press **[SET]** to confirm. If there is no operation, it will turn off automatically within 20 seconds
- After age setting, height sign will start flickering. Press **[▲]** or **[▼]** to adjust height, press **[SET]** to confirm. If there is no operation, it will turn off automatically within 20 seconds.

### 2 Start to weight your fat, water and muscle

- After all the data will be set, LCD shows "0.0", which means you can start to weigh your fat, water and muscle.
- Step on the scale, wait until the numbers stabilized, the final weight will be displayed and flash once.
- The result of the fat, water and muscles ranges will be displayed on the LCD one by one and these parameters will be showed three times in total and then automatically switched off.

## Change the unit

In the back of the scale, there is a "UNIT" button, press gently to change the unit. 3 units are selectable : KG, LB and ST.

**Notes:**

- . Gently tap the scale with your foot or hand, the LCD display shows the last record of the measurement from the user. (No record will appear in the first use).
- . The scale turns off automatically after 8 seconds.
- . Automatic shut off means that your personal parameters have been saved.

**Important notes:**

- While setting parameters you may set the scale to zero by tapping on the scale with your foot or your hand. Based on the current set or modified parameters you may begin measuring fat/water/muscles at any time.
- After completing the settings you may continue by pressing the SET button again and carry out another setup of all the parameters at any time.

**IMPORTANT SAFETY NOTICE**

- 1 Carefully unpack the product and take care not to throw away any part of the packaging material before you have found all parts of the product.
- 2 Under no circumstances should you repair or alter the product by yourself. Entrust all repairs and adjustments to a specialized company/service. Failure to do so may void the warranty.
- 3 Never use this product if you have wet feet. Danger of injuries from slipping.
- 4 Store and use the product away from the flammable or volatile materials and substances.
- 5 Don not expose the product to extreme heat, direct sunlight or excessive humidity and don not place in an overly dusty environment.
- 6 Do not place the product near heaters, open fires or other sources of heat.
- 7 This product is designed for household use. Don not use it in industrial surroundings of outdoors!
- 8 Do not use this product for other purposes than those intended.
- 9 Don not spray the product with water or any other liquid. Do not pour water or any other liquids into the product. Do not submerge the product in water or any other liquid.
- 10 Never jump on the scale.
- 11 Do not use rechargeable batteries in this product.
- 12 This scale is designed only for the weighing of people.
- 13 If you will not be using the product for an extended period of time, remove the batteries. Otherwise they may leak inside the product and cause irreparable damage.
- 14 The manufacturer takes no responsibility for damages caused by the incorrect use of this appliance and its accessories ( injuries, burns, scalding, fire, food spoilage, etc.).
- 15 Please do NOT scratch the paintings on the back of the glass which will destroy the art printing.

**CHILDREN AND FRAIL INDIVIDUALS**

- In order to ensure your children's safety, please keep all packaging (plastic bags, boxes, polystyrene etc.) out of their reach. **WARNING!** Do not allow small children to play with the foil as there is a danger of suffocation!
- This device is not intended to be used by individuals (including children) who have restricted physical, sensory or mental abilities and/or insufficient knowledge and/or experience, unless they are supervised by an individual who is responsible for their safety or have received instructions on how to use the device.
- Children should be supervised at all times in order to ensure that they do not play with the device.

## WARRANTY

1. This warranty is issued by EDC Poterek Sp. Jawna with address:  
ul. Poznańska 129/133, 05-850 Ożarów Mazowiecki, Phone: 022-7213571 w.10 – called hereinafter Esperanza.
2. This warranty confers upon you specific rights and remedies against Esperanza. Please read this document carefully, by making a claim under this warranty, you will be deemed to have understood and accept its terms.
3. Esperanza warrants that any Esperanza product together with attached accessories is free from defects in materials and workmanship under conditions of its normal use.
4. If you make a claim under this warranty, Esperanza will, at its option, repair or replace any defective product into new one by using new or use factory-refurbished parts. In case of using factory-refurbished parts they will be equal to new products.
5. Defective product should be returned together with receipt, commercial invoice or any other proof of purchase to the place of purchase. It should include all accessories attached in original packing. Esperanza has right to reject warranty if any accessories will be missed.
6. Esperanza is committed to respond to the complaint within 14 calendar days of its notification in the point of sale. Esperanza will repair the defective product in the shortest possible time. If will be need to import from abroad of any parts which are necessary to repair this time may be extended.
7. In case of replacing defective product onto new one, Esperanza has right to replace it with newer model with the same technical parameters if the old one model was replaced by new one during warranty period.
8. If you make claim under this warranty, you are deemed to have agreed that repair and (upon availability) replacement, as applicable, under the warranty services described herein, is your sole and exclusive remedy vis-a-vis Esperanza.

## EXCLUSIONS:

### This warranty does not extend to:

1. Any product not manufactured by Esperanza
2. Any product that has been damaged or rendered defective as a result of:
  - use it in other way than normal intended use
  - modification or internal interferes or changes of product
  - service or repair by anyone other than Esperanza service
  - improper store, transportation or packing
  - improper installation of product or software
  - mechanical, physical, chemical etc. damages
3. In no event shall Esperanza be liable for any damages, loss of business, profits, data or loss based on use Esperanza products with other devices.
4. Products with removed warranty and other labels and products spent or have hard use marks.

All products are covered by Esperanza's 2-year manufacturer's warranty while product is purchased by the consumer (the consumer is a person who buys product for the purpose of not professional or commercial activity). In the case of purchase for the company (for purposes relating to the activities of professional or commercial) warranty is 1 year from date of purchase. Esperanza products have all the necessary safety certificates required in the EU.

According to Directive 2012/19/EU, disposal of waste electrical and electronic equipment together with household waste is prohibited. Their components need to be recycled or disposed separately. Otherwise, they can cause the release of hazardous and noxious substances which pose a threat to the environment and human health. The law requires each consumer to return unnecessary or waste electrical and electronic equipment to the producers, sales agents or waste treatment facilities. Detailed conditions are governed by national legislation. The above symbol placed on the product, in the user manual or on the packaging refers to this requirement. Thanks to proper segregation and compliance with the provisions relating to waste disposal each user contributes significantly to environmental protection.

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